

photobiologyFilters Version 0.1.9

Catalogue of filters

Pedro J. Aphalo

July 5, 2014

1 Introduction

```
library(ggplot2)
library(photobiologyFilters)

## Loading required package: photobiology
## Loading required package: lubridate
## Loading required package: data.table
##
## Attaching package: 'data.table'
##
## The following objects are masked from 'package:lubridate':
##
##   hour, mday, month, quarter, wday, week, yday, year

library(photobiologygg)


## Loading required package: proto
## Loading required package: spls2R
## Loading required package: plyr
##
## Attaching package: 'plyr'
##
## The following object is masked from 'package:lubridate':
##
##   here

filter.plotter <- function(filter_name, w.low=280, w.high=800,
                             ylab="Spectral transmittance (%)"){
  spectrum.data <- data.frame(w.length=seq(280, 800, length.out=300))
  spectrum.data$transmittance <- calc_filter_multipliers(spectrum.data$w.length, filter_name, pc.out=TRUE)
  fig_linear <- ggplot(aes(x=w.length, y=transmittance), data=spectrum.data) +
    labs(x="Wavelength (nm)", y=ylab, title=filter_name) + ylim(0,100) +
    geom_line()
  # fig_log <- fig_linear + scale_y_log10(limits=c(1e-5,30))
  print(fig_linear + stat_peaks(span=71, ignore_threshold=0.25, colour="red", hjust=0, angle=90) +
        stat_valleys(span=51, ignore_threshold=-0.5, colour="blue", hjust=0, angle=90))
  # print(fig_log)
}
```

2 Dummy filters

2.1 Perfectly clear filter


```
filter.plotter("clear")
```





3 Plastic films


3.1 Cellulose diacetate




```
for (filter in c("acetate.115um.new", "acetate.250um.new", "acetate.480um.new")) {  
  filter.plotter(filter)  
}
```








```
for (filter in c("acetate.ageing0", "acetate.ageing20", "acetate.ageing30", "acetate.ageing60", "acetate.ageing100")) {  
  filter.plotter(filter)  
}
```









3.2 Polyester


```
filter.plotter("polyester.new")
```



3.3 Polythene



```
filter.plotter("polythene.new")  
filter.plotter("polythene.used")
```












3.4 Rosco theatrical filters

```
for (filter in c("clear.00.new", "uv.226.new", "uv.226.used", "canary.yellow.new",  
                "canary.yellow.used", "rose.pink.new", "neon.pink.used",  
                "moss.green.new", "moss.green.used")) {  
  filter.plotter(filter)  
}
```











3.5 Commercial greenhouse films from BPI Agri Visqueen



```
for (filter in c("solatrol.new", "luminance.new")) {  
  filter.plotter(filter)  
}
```





4 Plastic sheets


4.1 Plexiglas




```
for (filter in c("PLX0A000_XT", "PLX0A570_GT", "PLX0F00_GT", "PLX0Z023_GT")) {  
  filter.plotter(filter)  
}
```





```
for (filter in c("PLX1C33_GT", "PLX2C04_GT", "PLX3C01_GT", "PLX5C01_GT")) {
  filter.plotter(filter)
}
```











5 Optical glass filters

5.1 Schott long-pass filters

```
for (filter in c("gg395", "gg400", "gg435", "gg455", "gg475", "gg495")) {  
  filter.plotter(filter)  
}
```







5.2 Schott band-pass filters

```
for (filter in c("ug1", "ug11")) {
  filter.plotter(filter)
}
```

